


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	IOWH-032
<b>Catalog Number:</b>	BCC3922
<b>Cas No.:</b>	1191252-49-9
<b>Formula:</b>	C <sub>22</sub> H <sub>15</sub> Br <sub>2</sub> N <sub>3</sub> O <sub>4</sub>
<b>M.Wt:</b>	545.18
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	3-(3,5-dibromo-4-oxocyclohexa-2,5-dien-1-ylidene)-N-[(4-phenoxyphenyl)methyl]-1,2,4-oxadiazole-5-carboxamide
<b>SMILES:</b>	<chem>C1=CC=C(C=C1)OC2=CC=C(C=C2)CNC(=O)C3=NC(=C4C=C(C(=O)C(=C4)Br)Br)NO3</chem>
<b>Solubility:</b>	DMSO : ≥ 100 mg/mL (183.43 mM) *"≥" means soluble, but saturation unknown.



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893