


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	Sorafenib
<b>Catalog Number:</b>	BCN2174
<b>Cas No.:</b>	284461-73-0
<b>Formula:</b>	C <sub>21</sub> H <sub>16</sub> ClF <sub>3</sub> N <sub>4</sub> O <sub>3</sub>
<b>M.Wt:</b>	464.83
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Powder
<b>Synonyms:</b>	Bay 43-9006
<b>Chemical Name:</b>	4-[4-[[4-chloro-3-(trifluoromethyl)phenyl]carbamoylamino]phenoxy]-N-methylpyridine-2-carboxamide
<b>SMILES:</b>	<chem>CNC(=O)C1=NC=CC(=C1)OC2=CC=C(C=C2)NC(=O)NC3=CC(=C(C=C3)Cl)C(F)(F)F</chem>
<b>Solubility:</b>	DMSO : ≥ 45 mg/mL (96.81 mM) *"≥" means soluble, but saturation unknown.



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893