


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	Fraxiresinol 1-O-glucoside
<b>Catalog Number:</b>	BCN4439
<b>Cas No.:</b>	89199-94-0
<b>Formula:</b>	C <sub>27</sub> H <sub>34</sub> O <sub>13</sub>
<b>M.Wt:</b>	566.6
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Powder
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	(2S,3R,4S,5S,6R)-2-[[[(3R,3aS,6S,6aR)-3-(4-hydroxy-3,5-dimethoxyphenyl)-6-(4-hydroxy-3-methoxyphenyl)-3,4,6,6a-tetrahydro-1H-furo[3,4-c]furan-3a-yl]oxy]-6-(hydroxymethyl)oxane-3,4,5-triol
<b>SMILES:</b>	<chem>COC1=CC(=CC(=C1O)OC)C2C3(COC(C3CO2)C4=CC(=C(C=C4)O)OC)OC5C(C(C(C(O5)CO)O)O)O</chem>
<b>Solubility:</b>	N/A



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893