


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	CDK9 inhibitor
<b>Catalog Number:</b>	BCC1465
<b>Cas No.:</b>	1415559-43-1
<b>Formula:</b>	C <sub>26</sub> H <sub>21</sub> N <sub>5</sub> O <sub>4</sub> S
<b>M.Wt:</b>	499.54
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	N-[5-[[[6-[3-(1,3-dioxoisindol-2-yl)phenyl]pyrimidin-4-yl]amino]-2-methylphenyl]methanesulfonamide
<b>SMILES:</b>	<chem>CC1=C(C=C(C=C1)NC2=NC=NC(=C2)C3=CC(=CC=C3)N4C(=O)C5=CC=CC=C5C4=O)NS(=O)(=O)C</chem>
<b>Solubility:</b>	DMSO : 25 mg/mL (50.05 mM; Need ultrasonic) H <sub>2</sub> O : < 0.1 mg/mL (insoluble)



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15,Building B2, Tianfu Life Science Park,88 South Keyuan Road, Hi-Tech Zone,Chengdu,Sichuan 610041,PRC  
Tech Support:service@biocrick.com Tel:+86-28-8543-3893 FAX:+86-28-8543-3893