


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	Fimasartan
<b>Catalog Number:</b>	BCC5552
<b>Cas No.:</b>	247257-48-3
<b>Formula:</b>	C <sub>27</sub> H <sub>31</sub> N <sub>7</sub> O <sub>5</sub>
<b>M.Wt:</b>	501.65
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	2-[2-butyl-4-methyl-6-oxo-1-[[4-[2-(2H-tetrazol-5-yl)phenyl]phenyl]methyl]pyrimidin-5-yl]-N,N-dimethylethanethioamide
<b>SMILES:</b>	<chem>CCCC1=NC(=C(C(=O)N1CC2=CC=C(C=C2)C3=CC=CC=C3C4=NNN=C4)CC(=S)N(C)C)C</chem>
<b>Solubility:</b>	DMSO : ≥ 49 mg/mL (97.68 mM) *"≥" means soluble, but saturation unknown.



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893