


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	BMS 193885
<b>Catalog Number:</b>	BCC7613
<b>Cas No.:</b>	679839-66-8
<b>Formula:</b>	C36H48N4O9
<b>M.Wt:</b>	680.79
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	dimethyl 4-[3-[3-[4-(3-methoxyphenyl)piperidin-1-yl]propylcarbamoylamino]phenyl]-2,6-dimethyl-1,4-dihydropyridine-3,5-dicarboxylate;(2S)-2-hydroxypropionic acid
<b>SMILES:</b>	<chem>CC1=C(C(C(=C(N1)C)C(=O)OC)C2=CC(=CC=C2)NC(=O)NCCCN3CCC(CC3)C4=CC(=CC=C4)OC)C(=O)OC.CC(C(=O)O)O</chem>
<b>Solubility:</b>	Soluble to 20 mM in water and to 100 mM in DMSO



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893