


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	3,19-Dihydroxy-6,23-dioxo-12-ursen-28-oic acid
<b>Catalog Number:</b>	BCN1471
<b>Cas No.:</b>	261768-88-1
<b>Formula:</b>	C <sub>30</sub> H <sub>44</sub> O <sub>6</sub>
<b>M.Wt:</b>	500.7
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Powder
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	(1R,2R,4aS,6aR,6aS,6bR,8aR,9S,10S,12aR,14bS)-9-formyl-1,10-dihydroxy-1,2,6a,6b,9,12a-hexamethyl-8-oxo-3,4,5,6,6a,7,8a,10,11,12,13,14b-dodecahydro-2H-picene-4a-carboxylic acid
<b>SMILES:</b>	<chem>CC1CCC2(CCC3(C(=CCC4C3(CC(=O)C5C4(CCC(C5(C)C=O)O)C)C)C2C1(C)O)C)C(=O)O</chem>
<b>Solubility:</b>	N/A



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15,Building B2, Tianfu Life Science Park,88 South Keyuan Road, Hi-Tech Zone,Chengdu,Sichuan 610041,PRC  
Tech Support:service@biocrick.com Tel:+86-28-8543-3893 FAX:+86-28-8543-3893