


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	3-O-cis-p-Coumaroyltormentic acid
<b>Catalog Number:</b>	BCN3184
<b>Cas No.:</b>	121072-40-0
<b>Formula:</b>	C39H54O7
<b>M.Wt:</b>	634.9
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Powder
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	(1R,2R,4aS,6aS,6bR,10R,11R,12aR,14bS)-1,11-dihydroxy-10-[(Z)-3-(4-hydroxyphenyl)prop-2-enoyl]oxy-1,2,6a,6b,9,9,12a-heptamethyl-2,3,4,5,6,6a,7,8,8a,10,11,12,13,14b-tetradecahydronicene-4a-carboxylic acid
<b>SMILES:</b>	<chem>CC1CCC2(CCC3(C(=CCC4C3(CCC5C4(CC(C(C5(C)C)OC(=O)C=CC6=C(C=C(C=C6)O)O)C)C)C2C1(C)O)C)C(=O)O</chem>
<b>Solubility:</b>	N/A



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893