


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	Orthosiphol A
<b>Catalog Number:</b>	BCN8879
<b>Cas No.:</b>	142741-25-1
<b>Formula:</b>	C38H44O11
<b>M.Wt:</b>	676.7
<b>Purity:</b>	>98%
<b>Physical Description:</b>	Powder
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	[(2R,4R,4aS,4bS,5R,6S,7S,8aS,10R,10aR)-6,10-diacetyloxy-5-benzoyloxy-2-ethenyl-7,10a-dihydroxy-2,4b,8,8-tetramethyl-1-oxo-4,4a,5,6,7,8a,9,10-octahydro-3H-phenanthren-4-yl] benzoate
<b>SMILES:</b>	<chem>CC(=O)OC1CC2C(C(C(C(C2(C3C1(C(=O)C(CC3OC(=O)C4=CC=CC=C4)(C)C=C)O)C)OC(=O)C5=CC=CC=C5)OC(=O)C)O)(C)C</chem>
<b>Solubility:</b>	N/A



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: [service@biocrick.com](mailto:service@biocrick.com) Tel: +86-28-8543-3893 FAX: +86-28-8543-3893