

### [ Analysis Information ]

Product Name:	COG 133
Catalog No.:	BCC1047
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:	Venusil ASB C18,4.6*250mm,5 $\mu$ m,300A	
Solvent A:	0.1%Trifluoroacetic in 100% Acetonitrile	
Solvent B:	0.1%Trifluoroacetic in 100% Water	
Gradient:	Time	A
	0.0min	10%
	30.0min	70%
	35.0min	100%
	35.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	220nm	
Volume:	15 $\mu$ l	

### [ Result ]

