

### [ Analysis Information ]

|               |              |
|---------------|--------------|
| Product Name: | 10Panx       |
| Catalog No.:  | BCC1245      |
| Operator:     | Jianhong Lee |
| Batch No.:    | b01          |

### [ Analysis Condition ]

Column: Venusil ASB C18, 4.6\*250mm, 5μm, 300A

Solvent A: 0.1% Trifluoroacetic in 100% Acetonitrile

Solvent B: 0.1% Trifluoroacetic in 100% Water

| Gradient: | Time    | A    |
|-----------|---------|------|
|           | 0.0min  | 10%  |
|           | 30.0min | 70%  |
|           | 35.0min | 100% |
|           | 35.0min | Stop |

Flow rate: 0.5ml/min

Wavelength: 220nm

Volume: 15μl

### [ Result ]

