

### [ Analysis Information ]

|               |                  |
|---------------|------------------|
| Product Name: | VS-5584 (SB2343) |
| Catalog No.:  | BCC2047          |
| Operator:     | Jianhong Lee     |
| Batch No.:    | b01              |

### [ Analysis Condition ]

|             |   |      |
|-------------|---|------|
| Column:     | Waters XBridge C18, 3 $\mu$ m, 2.1 mm $\times$ 100 mm |      |
| Solvent A:  | 0.1% H3PO4 in 100% Acetonitrile                       |      |
| Solvent B:  | 0.1% H3PO4 in 100% Water                              |      |
| Gradient:   | Time  | A    |
|             | 0.0min  | 10%  |
|             | 3.5min  | 95%  |
|             | 8.0min  | 95%  |
|             | 9.0min  | 10%  |
|             | 10.0min   | 10%  |
|             | 10.0min   | Stop |
| Flow rate:  | 0.5ml/min   |      |
| Wavelength: | 254nm   |      |
| Volume:     | 15 $\mu$ l  |      |

### [ Result ]

