


## PHYSICAL AND CHEMICAL PROPERTIES

|                              |  |
|------------------------------|--|
| <b>Product Name:</b>         | XMD17-109  |
| <b>Catalog Number:</b>       | BCC2061  |
| <b>Cas No.:</b>              | 1435488-37-1   |
| <b>Formula:</b>              | C <sub>36</sub> H <sub>46</sub> N <sub>8</sub> O <sub>3</sub>  |
| <b>M.Wt:</b>                 | 638.8  |
| <b>Purity:</b>               | >98%   |
| <b>Physical Description:</b> | N/A  |
| <b>Synonyms:</b>             | N/A  |
| <b>Chemical Name:</b>        | 11-cyclopentyl-2-[2-ethoxy-4-[4-(4-methylpiperazin-1-yl)piperidine-1-carbonyl]anilino]-5-methylpyrimido[4,5-b][1,4]benzodiazepin-6-one |
| <b>SMILES:</b>               | <chem>CCOC1=C(C=CC(=C1)C(=O)N2CCC(CC2)N3CCN(CC3)C)NC4=NC=C5C(=N4)N(C6=CC=CC=C6C(=O)N5C)C7CCCC7</chem>                                  |
| <b>Solubility:</b>           | DMSO : ≥ 100 mg/mL (156.54 mM)<br>*"≥" means soluble, but saturation unknown.  |



bioctick.com

## STORAGE AND SHIPPING

|                            |   |
|----------------------------|---|
| <b>Storage:</b>            | Store the product in sealed, cool and dry condition.  |
| <b>General tips:</b>       | For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months. |
| <b>Shipping Condition:</b> | Ship with RT , or blue ice upon request   |

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893