

[Analysis Information]

Product Name:	Aspirin (Acetylsalicylic acid)
Catalog No.:	BCC2097
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18, 3 μ m, 2.1mm \times 100 mm										
Solvent A:	0.1% H ₃ PO ₄ in 100% Acetonitrile										
Solvent B:	0.1% H ₃ PO ₄ in 100% Water										
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>10%</td></tr><tr><td>6.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	10%	6.0min	95%	10.0min	95%	10.0min	Stop
Time	A										
0.0min	10%										
6.0min	95%										
10.0min	95%										
10.0min	Stop										
Flow rate:	0.5ml/min										
Wavelength:	230nm										
Volume:	15 μ l										

[Result]

