

[Analysis Information]

Product Name: Danoprevir (RG7227)
Catalog No.: BCC2106
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

| Time | A |
|--------|------|
| 0.0min | 10% |
| 6.0min | 95% |
| 10.0 | 95% |
| 10.0 | Stop |

Flow rate: 0.5ml/min

Wavelength: 210nm

Volume: 15 μ l

[Result]

