

[Analysis Information]

Product Name:	LAQ824 (NVP-LAQ824,Dacinostat)
Catalog No.:	BCC2160
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18, 3 μ m, 2.1mm \times 100 mm	
Solvent A:	0.1% H3PO4 in 100% Acetonitrile	
Solvent B:	0.1% H3PO4 in 100% Water	
Gradient:	Time	A
	0.0min	0%
	2.0min	0%
	6.0min	95%
	9.0min	95%
	10.0min	0%
	10.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	254nm	
Volume:	15 μ l	

[Result]

