

### [ Analysis Information ]

Product Name:	Scriptaid
Catalog No.:	BCC2163
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column: Waters XBridge C18, 3  $\mu$ m, 3 mm  $\times$  100 mm

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

Time	A
0.0min	95%
5.0min	95%
10.0min	95%
10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 254nm

Volume: 15 $\mu$ l

### [ Result ]

