

[Analysis Information]

Product Name:	Barasertib (AZD1152-HQPA)
Catalog No.:	BCC2168
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18 (3 μ m, 2.1 mm \times 100 mm)												
Solvent A:	0.1% H ₃ PO ₄ in 100% Acetonitrile												
Solvent B:	0.1% H ₃ PO ₄ in 100% Water												
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>5%</td></tr><tr><td>8.0min</td><td>95%</td></tr><tr><td>9.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	5%	8.0min	95%	9.0min	95%	10.0min	95%	10.0min	Stop
Time	A												
0.0min	5%												
8.0min	95%												
9.0min	95%												
10.0min	95%												
10.0min	Stop												
Flow rate:	0.5ml/min												
Wavelength:	254nm												
Volume:	15 μ l												

[Result]

