

### [ Analysis Information ]

Product Name:	Tofacitinib (CP-690550) Citrate
Catalog No.:	BCC2189
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:	Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm														
Solvent A:	0.1% H3PO4 in 100% Acetonitrile														
Solvent B:	0.1% H3PO4 in 100% Water														
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>0%</td></tr><tr><td>2.0min</td><td>0%</td></tr><tr><td>6.0min</td><td>95%</td></tr><tr><td>9.0min</td><td>95%</td></tr><tr><td>10.0min</td><td>0%</td></tr><tr><td>10.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	0%	2.0min	0%	6.0min	95%	9.0min	95%	10.0min	0%	10.0min	Stop
Time	A														
0.0min	0%														
2.0min	0%														
6.0min	95%														
9.0min	95%														
10.0min	0%														
10.0min	Stop														
Flow rate:	0.5ml/min														
Wavelength:	210nm														
Volume:	15 $\mu$ l														

### [ Result ]

