

[Analysis Information]

Product Name:	Tenovin-1
Catalog No.:	BCC2239
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column: Waters XBridge C18, 3 μ m, 2.1 mm \times 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

Time	A
0.0min	10%
3.5min	95%
8.0min	95%
9.0min	10%
10.0min	10%
10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 254nm

Volume: 15 μ l

[Result]

