

[Analysis Information]

Product Name: WY-14643 (Pirinixic Acid)
Catalog No.: BCC2265
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm
Solvent A: 0.1% H3PO4 in 100% Acetonitrile
Solvent B: 0.1% H3PO4 in 100% Water
Gradient:

Time	A
0.0min	10%
3.5min	95%
8.0min	95%
9.0min	10%
10.0min	10%
10.0min	Stop

Flow rate: 0.5ml/min
Wavelength: 270nm
Volume: 15 μ l

[Result]

