

### [ Analysis Information ]

|               |              |
|---------------|--------------|
| Product Name: | Repaglinide  |
| Catalog No.:  | BCC2504      |
| Operator:     | Jianhong Lee |
| Batch No.:    | b01          |

### [ Analysis Condition ]

|             |  |      |
|-------------|--|------|
| Column:     | Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm |      |
| Solvent A:  | 0.1% H3PO4 in 100% Acetonitrile              |      |
| Solvent B:  | 0.1% H3PO4 in 100% Water                     |      |
| Gradient:   | Time   | A    |
|             | 0.0min                                       | 0%   |
|             | 2.0min                                       | 0%   |
|             | 6.0min                                       | 95%  |
|             | 9.0min                                       | 95%  |
|             | 10.0min                                      | 0%   |
|             | 10.0min                                      | Stop |
| Flow rate:  | 0.5ml/min                                    |      |
| Wavelength: | 254nm  |      |
| Volume:     | 15 $\mu$ l                                   |      |

### [ Result ]

