

[Analysis Information]

Product Name:	Fmoc-L-Arg(Aloc)2-OH
Catalog No.:	BCC2564
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:					
Solvent A:					
Solvent B:					
Gradient:	<table border="1"> <tr> <th>Time</th> <th>A</th> </tr> <tr> <td></td> <td>Stop</td> </tr> </table>	Time	A		Stop
Time	A				
	Stop				
Flow rate:	0.5ml/min				
Wavelength:					
Volume:	15µl				

[Result]

