

[Analysis Information]

| | |
|---------------|--------------|
| Product Name: | Fmoc-Ser-OH |
| Catalog No.: | BCC3541 |
| Operator: | Jianhong Lee |
| Batch No.: | b01 |

[Analysis Condition]

| | | |
|-------------|--|------|
| Column: | Athena C18, 3 μ m, 2.1mm \times 100 mm | |
| Solvent A: | 0.1% H ₃ PO ₄ in 100% Acetonitrile | |
| Solvent B: | 0.1% H ₃ PO ₄ in 100% Water | |
| Gradient: | Time | A |
| | 0.0min | 10% |
| | 6.0min | 95% |
| | 10.0min | 95% |
| | 10.0min | Stop |
| Flow rate: | 0.5ml/min | |
| Wavelength: | 210nm | |
| Volume: | 15 μ l | |

[Result]

