

[Analysis Information]

Product Name:	Panobinostat (LBH589)
Catalog No.:	BCC3601
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18 (3 μ m, 2.1 mm \times 100 mm)	
Solvent A:	0.1% H ₃ PO ₄ in 100% Acetonitrile	
Solvent B:	0.1% H ₃ PO ₄ in 100% Water	
Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	254nm	
Volume:	15 μ l	

[Result]

