

[Analysis Information]

Product Name:	PSI-6206
Catalog No.:	BCC3609
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Waters XBridge C18, 3 μ m, 3 mm \times 100 mm	
Solvent A:	0.1% H3PO4 in 100% Acetonitrile	
Solvent B:	0.1% H3PO4 in 100% Water	
Gradient:	Time	A
	0.0min	0%
	3.0min	0%
	7.0min	35%
	12.0min	95%
	12.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	254nm	
Volume:	15 μ l	

[Result]

