

[Analysis Information]

Product Name:	Semagacestat (LY450139)
Catalog No.:	BCC3610
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	xBridge C18,4.6*150mm,3.5 μ m								
Solvent A:	0.1%Trifluoroacetic in 100% Acetonitrile								
Solvent B:	0.1%Trifluoroacetic in 100% Water 0.1%Trifluoroacetic in 100% Water								
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>20%</td></tr><tr><td>16.0min</td><td>80%</td></tr><tr><td>16.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	20%	16.0min	80%	16.0min	Stop
Time	A								
0.0min	20%								
16.0min	80%								
16.0min	Stop								
Flow rate:	0.5ml/min								
Wavelength:	210nm								
Volume:	15 μ l								

[Result]

