

[Analysis Information]

Product Name: Lapatinib
Catalog No.: BCC3633
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm
Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile
Solvent B: 0.1% H₃PO₄ in 100% Water

| Gradient: | Time | A |
|-----------|---------|------|
| | 0.0min | 0% |
| | 2.0min | 0% |
| | 7.0min | 95% |
| | 10.0min | 95% |
| | 10.0min | Stop |

Flow rate: 0.5ml/min
Wavelength: 210nm
Volume: 15 μ l

[Result]

