

### [ Analysis Information ]

Product Name:	BX795
Catalog No.:	BCC3635
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column: Athena C18 (3  $\mu$ m, 2.1 mm  $\times$  100 mm)

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 220nm

Volume: 15 $\mu$ l

### [ Result ]

