

### [ Analysis Information ]

Product Name: StemRegenin 1 (SR1)  
Catalog No.: BCC3637  
Operator: Jianhong Lee  
Batch No.: b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm  
Solvent A: 0.1% H3PO4 in 100% Acetonitrile  
Solvent B: 0.1% H3PO4 in 100% Water  
Gradient:

Time	A
0.0min	0%
2.0min	0%
7.0min	95%
10.0min	95%
10.0min	Stop

Flow rate: 0.5ml/min  
Wavelength: 319nm  
Volume: 15 $\mu$ l

### [ Result ]

