

### [ Analysis Information ]

Product Name: Regorafenib  
Catalog No.: BCC3646  
Operator: Jianhong Lee  
Batch No.: b01

### [ Analysis Condition ]

Column: Waters XBridge C18, 3  $\mu$ m, 2.1 mm  $\times$  100 mm  
Solvent A: 0.1% H3PO4 in 100% Acetonitrile  
Solvent B: 0.1% H3PO4 in 100% Water

Time	A
0.0min	10%
3.5min	95%
8.0min	95%
9.0min	10%
10.0min	10%
10.0min	Stop

Flow rate: 0.5ml/min  
Wavelength: 220nm  
Volume: 15 $\mu$ l

### [ Result ]

