

### [ Analysis Information ]

Product Name:	TAE684 (NVP-TAE684)
Catalog No.:	BCC3660
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:	Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm	
Solvent A:	0.1% H3PO4 in 100% Acetonitrile	
Solvent B:	0.1% H3PO4 in 100% Water	
Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	254nm	
Volume:	15 $\mu$ l	

### [ Result ]

