

[Analysis Information]

| | |
|---------------|--------------|
| Product Name: | Deferasirox |
| Catalog No.: | BCC3924 |
| Operator: | Jianhong Lee |
| Batch No.: | b01 |

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm

Solvent A: 0.1% H₃PO₄ in 100% Acetonitrile

Solvent B: 0.1% H₃PO₄ in 100% Water

| Time | A |
|---------|------|
| 0.0min | 10% |
| 3.5min | 95% |
| 8.0min | 95% |
| 9.0min | 10% |
| 10.0min | 10% |
| 10.0min | Stop |

Flow rate: 0.5ml/min

Wavelength: 299nm

Volume: 15 μ l

[Result]

