

[Analysis Information]

Product Name:	ZIP
Catalog No.:	BCC4003
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Venusil ASB C18,4.6*250mm,5µm,300A 0.1%Trifluoroacetic in 100% Acetonitrile	
Solvent A:	0.1%Trifluoroacetic in 100% Water	
Solvent B:		
Gradient:	Time	A
	0.0min	10%
	20.0min	70%
	35.0min	100%
	35.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	220nm	
Volume:	15µl	

[Result]

