

[Analysis Information]

Product Name: ZIP
 Catalog No.: BCC4003
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: Venusil ASB C18, 4.6*250mm, 5 μ m, 300A 0.1% Trifluoroacetic in 100% Acetonitrile
 Solvent A: 0.1% Trifluoroacetic in 100% Water
 Solvent B:
 Gradient:

Time	A
0.0min	10%
20.0min	70%
35.0min	100%
35.0min	Stop

Flow rate: 0.5ml/min
 Wavelength: 220nm
 Volume: 15 μ l

[Result]

