



## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	AST-1306 TsOH
<b>Catalog Number:</b>	BCC4043
<b>Cas No.:</b>	1050500-29-2
<b>Formula:</b>	C <sub>31</sub> H <sub>26</sub> ClFN <sub>4</sub> O <sub>5</sub> S
<b>M.Wt:</b>	621.08
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	AST-1306 (p-Toluenesulfonic acid); AST1306 (p-Toluenesulfonic acid); AST 1306 (p-Toluenesulfonic acid)
<b>Chemical Name:</b>	N-[4-[3-chloro-4-[(3-fluorophenyl)methoxy]anilino]quinazolin-6-yl]prop-2-enamide; 4-methylbenzenesulfonic acid
<b>SMILES:</b>	<chem>CC1=CC=C(C=C1)S(=O)(=O)O.C=CC(=O)NC1=CC2=C(C=C1)N=CN=C2NC3=CC(=C(C=C3)OCC4=CC(=CC=C4)F)Cl</chem>
<b>Solubility:</b>	DMSO : 50 mg/mL (80.50 mM; Need ultrasonic) H <sub>2</sub> O : < 0.1 mg/mL (insoluble)



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure*

*that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893