

### [ Analysis Information ]

Product Name:	Doripenem
Catalog No.:	BCC4094
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column: Athena C18, 3  $\mu$ m, 2.1mm  $\times$  100 mm

Solvent A: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Acetonitrile

Solvent B: 0.1% H<sub>3</sub>PO<sub>4</sub> in 100% Water

Gradient:	Time	A
	0.0min	0%
	2.0min	0%
	7.0min	95%
	10.0	95%
	10.0	Stop

Flow rate: 0.5ml/min

Wavelength: 310nm

Volume: 15 $\mu$ l

### [ Result ]

