

[Analysis Information]

Product Name: CC-401
Catalog No.: BCC4269
Operator: Jianhong Lee
Batch No.: b01

[Analysis Condition]

Column: Athena C18, 3 μ m, 2.1mm \times 100 mm
Solvent A: 0.1% H3PO4 in 100% Acetonitrile
Solvent B: 0.1% H3PO4 in 100% Water
Gradient:

| Time | A |
|---------|------|
| 0.0min | 0% |
| 2.0min | 0% |
| 6.0min | 95% |
| 9.0min | 95% |
| 10.0min | 0% |
| 10.0min | Stop |

Flow rate: 0.5ml/min
Wavelength: 254nm
Volume: 15 μ l

[Result]

