

### [ Analysis Information ]

Product Name:	Irinotecan HCl Trihydrate
Catalog No.:	BCC5091
Operator:	Jianhong Lee
Batch No.:	b01

### [ Analysis Condition ]

Column:	Athena C18, 3 $\mu$ m, 2.1mm $\times$ 100 mm										
Solvent A:	0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Acetonitrile										
Solvent B:	0.1% H <sub>3</sub> PO <sub>4</sub> in 100% Water										
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>5%</td></tr><tr><td>9.0min</td><td>95%</td></tr><tr><td>15.0min</td><td>95%</td></tr><tr><td>15.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	5%	9.0min	95%	15.0min	95%	15.0min	Stop
Time	A										
0.0min	5%										
9.0min	95%										
15.0min	95%										
15.0min	Stop										
Flow rate:	0.5ml/min										
Wavelength:	254nm										
Volume:	15 $\mu$ l										

### [ Result ]

