


## PHYSICAL AND CHEMICAL PROPERTIES

<b>Product Name:</b>	PPNDS
<b>Catalog Number:</b>	BCC7015
<b>Cas No.:</b>	1021868-77-8
<b>Formula:</b>	C <sub>18</sub> H <sub>11</sub> N <sub>4</sub> Na <sub>4</sub> O <sub>14</sub> PS <sub>2</sub>
<b>M.Wt:</b>	694.36
<b>Purity:</b>	>98%
<b>Physical Description:</b>	N/A
<b>Synonyms:</b>	N/A
<b>Chemical Name:</b>	tetrasodium;3-[[4-formyl-5-hydroxy-6-methyl-3-(phosphonatooxymethyl)pyridin-2-yl]diazenyl]-7-nitronaphthalene-1,5-disulfonate
<b>SMILES:</b>	<chem>CC1=C(C(=C(C(=N1)N=NC2=CC(=C3C=C(C=C(C3=C2)S(=O)(=O)[O-]))N+)(=O)[O-])S(=O)(=O)[O-])COP(=O)([O-])[O-])C=O)O.[Na+].[Na+].[Na+].[Na+]</chem>
<b>Solubility:</b>	Soluble to 100 mM in water



## STORAGE AND SHIPPING

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15,Building B2, Tianfu Life Science Park,88 South Keyuan Road, Hi-Tech Zone,Chengdu,Sichuan 610041,PRC  
Tech Support:service@biocrick.com Tel:+86-28-8543-3893 FAX:+86-28-8543-3893