

[Analysis Information]

Product Name:	Gefitinib
Catalog No.:	BCN2173
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column: Venusil XBridge C18 ,4.6*150mm, 3.5um,300A

Solvent A: 0.1%Trifluoroacetic in 100% Acetonitrile

Solvent B: 0.1%Trifluoroacetic in 100% Water

Time	A
0.0min	95%
3.0min	95%
6.0min	95%
6.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 262nm

Volume: 15µl

[Result]

