

[Analysis Information]

Product Name:	Caffeic acid phenethyl ester
Catalog No.:	BCN2695
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18, 3 μ m, 2.1mm \times 100 mm														
Solvent A:	0.1% H3PO4 in 100% Acetonitrile														
Solvent B:	0.1% H3PO4 in 100% Water														
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>10%</td></tr><tr><td>3.5min</td><td>95%</td></tr><tr><td>8.0min</td><td>95%</td></tr><tr><td>9.0min</td><td>10%</td></tr><tr><td>10.0min</td><td>10%</td></tr><tr><td>10.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	10%	3.5min	95%	8.0min	95%	9.0min	10%	10.0min	10%	10.0min	Stop
Time	A														
0.0min	10%														
3.5min	95%														
8.0min	95%														
9.0min	10%														
10.0min	10%														
10.0min	Stop														
Flow rate:	0.5ml/min														
Wavelength:	250nm														
Volume:	15 μ l														

[Result]

