

## PHYSICAL AND CHEMICAL PROPERTIES

**Product Name:** Picfeltarraenin X

**Catalog Number:** BCN2859

**Cas No.:** 1391826-61-1

**Formula:** C<sub>36</sub>H<sub>54</sub>O<sub>11</sub>

**M.Wt:** 662.81

**Purity:** >98%


**Physical Description:** Powder

**Synonyms:** N/A

**Chemical Name:** (2S)-2-[(2S,3S,8S,9R,10R,13R,14S,16R,17R)-3,16-dihydroxy-4,4,9,13,14-pentamethyl-11-oxo-2-[(2R,3R,4S,5S,6R)-3,4,5-trihydroxy-6-(hydroxymethyl)oxan-2-yl]oxy-1,2,3,7,8,10,12,15,16,17-decahydrocyclopenta[a]phenanthren-17-yl]-2-methyl-5-propan-2-ylfuran-3-one

**SMILES:** CC(C)C1=CC(=O)C(O1)(C)C2C(CC3(C2(CC(=O)C4(C3CC=C5C4CC(C(C5(C)C)O)OC6C(C(C(C(O6)CO)O)O)O)C)C)C)O

**Solubility:** N/A



## STORAGE AND SHIPPING

**Storage:** Store the product in sealed, cool and dry condition.

**General tips:** For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.

**Shipping Condition:** Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure*

*that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893