


## PHYSICAL AND CHEMICAL PROPERTIES

|                              |   |
|------------------------------|---|
| <b>Product Name:</b>         | Myricetin 3-O-beta-D-xylopyranosyl(1-2)-beta-D-glucopyranoside  |
| <b>Catalog Number:</b>       | BCN8140   |
| <b>Cas No.:</b>              | 142449-93-2   |
| <b>Formula:</b>              | C <sub>26</sub> H <sub>28</sub> O <sub>17</sub>   |
| <b>M.Wt:</b>                 | 612.5   |
| <b>Purity:</b>               | >98%  |
| <b>Physical Description:</b> | Powder  |
| <b>Synonyms:</b>             | N/A   |
| <b>Chemical Name:</b>        | 3-[(2S,3R,4S,5S,6R)-4,5-dihydroxy-6-(hydroxymethyl)-3-[(2S,3R,4S,5R)-3,4,5-trihydroxyoxan-2-yl]oxyoxan-2-yl]oxy-5,7-dihydroxy-2-(3,4,5-trihydroxyphenyl)chromen-4-one |
| <b>SMILES:</b>               | <chem>C1C(C(C(C(O1)OC2C(C(C(OC2OC3=C(OC4=CC(=CC(=C4C3=O)O)O)O)C5=CC(=C(C(=C5)O)O)O)CO)O)O)O)O</chem>  |
| <b>Solubility:</b>           | N/A   |



## STORAGE AND SHIPPING

|                            |   |
|----------------------------|---|
| <b>Storage:</b>            | Store the product in sealed, cool and dry condition.  |
| <b>General tips:</b>       | For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months. |
| <b>Shipping Condition:</b> | Ship with RT , or blue ice upon request   |

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893