

[Analysis Information]

Product Name:	Acetyl Angiotensinogen (1-14), porcine
Catalog No.:	BCC1023
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Venusil ASB C18,4.6*250mm,5µm,300A										
Solvent A:	0.1%Trifluoroacetic in 100% Acetonitrile										
Solvent B:	0.1%Trifluoroacetic in 100% Water										
Gradient:	<table border="1"><thead><tr><th>Time</th><th>A</th></tr></thead><tbody><tr><td>0.0min</td><td>10%</td></tr><tr><td>30.0min</td><td>70%</td></tr><tr><td>35.0min</td><td>100%</td></tr><tr><td>35.0min</td><td>Stop</td></tr></tbody></table>	Time	A	0.0min	10%	30.0min	70%	35.0min	100%	35.0min	Stop
Time	A										
0.0min	10%										
30.0min	70%										
35.0min	100%										
35.0min	Stop										
Flow rate:	0.5ml/min										
Wavelength:	220nm										
Volume:	15µl										

[Result]

