

[Analysis Information]

| | |
|---------------|------------------------|
| Product Name: | Enzastaurin (LY317615) |
| Catalog No.: | BCC1100 |
| Operator: | Jianhong Lee |
| Batch No.: | b01 |

[Analysis Condition]

| | | |
|-------------|--|------|
| Column: | Athena C18 (3 μ m, 2.1 mm \times 100 mm) | |
| Solvent A: | 0.1% H ₃ PO ₄ in 100% Acetonitrile | |
| Solvent B: | 0.1% H ₃ PO ₄ in 100% Water | |
| Gradient: | Time | A |
| | 0.0min | 10% |
| | 3.5min | 95% |
| | 8.0min | 95% |
| | 9.0min | 10% |
| | 10.0min | 10% |
| | 10.0min | Stop |
| Flow rate: | 0.5ml/min | |
| Wavelength: | 254nm | |
| Volume: | 15 μ l | |

[Result]

