

[Analysis Information]

Product Name: Q-VD-OPh hydrate
 Catalog No.: BCC1125
 Operator: Jianhong Lee
 Batch No.: b01

[Analysis Condition]

Column: xBridge C18,4.6*150mm,3.5μ m
 Solvent A: 0.1%Trifluoroacetic in 100% Acetonitrile
 Solvent B: 0.1%Trifluoroacetic in 100% Water
 Gradient:

Time	A
0.0min	30%
A	15.0min
30%	70%
15.0min	B
70%	70%
B	30%
70%	0.6ml/min
30%	254nm
	10μl

Rank
 RetTime(min)
 Area(mAU's)

Area
 %
 5.941
 8.074
 8.306
 9.820
 1030074
 13338587
 12050175
 231309
 3.87
 50.05
 45.22
 0.87
 Stop

Flow rate: 0.5ml/min

Manufactured By BioCrick Biotech - Not Fully Tested - Research Use Only - Not For Human or Veterinary Use

Volume: 15μl

Email: info@biocrick.com