

[Analysis Information]

Product Name:	MEK162 (ARRY-162, ARRY-438162)
Catalog No.:	BCC1148
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column:	Athena C18, 3 μ m, 2.1mm $\times$ 100 mm	
Solvent A:	0.1% H3PO4 in 100% Acetonitrile	
Solvent B:	0.1% H3PO4 in 100% Water	
Gradient:	Time	A
	0.0min	10%
	3.5min	95%
	8.0min	95%
	9.0min	10%
	10.0min	10%
	10.0min	Stop
Flow rate:	0.5ml/min	
Wavelength:	258nm	
Volume:	15 μ l	

[Result]

