


---

## PHYSICAL AND CHEMICAL PROPERTIES

---

<b>Product Name:</b>	NSC 23766	
<b>Catalog Number:</b>	BCC1149	
<b>Cas No.:</b>	1177865-17-6	
<b>Formula:</b>	C <sub>24</sub> H <sub>38</sub> Cl <sub>3</sub> N <sub>7</sub>	
<b>M.Wt:</b>	530.96	
<b>Purity:</b>	>98%	
<b>Physical Description:</b>	N/A	
<b>Synonyms:</b>	N/A	
<b>Chemical Name:</b>	6-N-[2-[5-(diethylamino)pentan-2-ylamino]-6-methylpyrimidin-4-yl]-2-methylquinoline-4,6-diamine;trihydrochloride	
<b>SMILES:</b>	CCN(CC)CCCC(C)NC1=NC(=CC(=N1)NC2=CC3=C(C=C2)N=C(C=C3N)C)C.Cl.Cl.Cl	
<b>Solubility:</b>	DMSO : 33.33 mg/mL (62.77 mM; Need ultrasonic) H <sub>2</sub> O : ≥ 32 mg/mL (60.27 mM) *"≥" means soluble, but saturation unknown.	

---

## STORAGE AND SHIPPING

---

<b>Storage:</b>	Store the product in sealed, cool and dry condition.
<b>General tips:</b>	For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.
<b>Shipping Condition:</b>	Ship with RT , or blue ice upon request

---

## INTENDED USE

---

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

---

## CAUTION!

---

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure*

*that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893