

[Analysis Information]

Product Name:	MG-115
Catalog No.:	BCC1237
Operator:	Jianhong Lee
Batch No.:	b01

[Analysis Condition]

Column: Venusil ASB C18, 4.6*250mm, 5μm, 300A

Solvent A: 0.1% Trifluoroacetic in 100% Acetonitrile

Solvent B: 0.1% Trifluoroacetic in 100% Water

Gradient:	Time	A
	0.0min	10%
	30.0min	70%
	35.0min	100%
	35.0min	Stop

Flow rate: 0.5ml/min

Wavelength: 220nm

Volume: 15μl

[Result]

